

making pH measurement SMART



no stand needed neasure

- Familiar Android tablet platform
- Live display of single and multiple pH electrodes with one tablet
- Uses any standard pH electrode with an S7 Connector
- Visible calibration information
- Easy data storage and sharing
- Warning alerts and alarms
- Stores up to 10,000 readings
- Automated fault diagnosis
- Electrode health reminder

anywhere



patented system







6 meters in 1 Choose between one large single screen meter or simultaneously view the readings of up to 6 SMART Caps at once on a single tablet.



The TRUEpH SMART Cap from TRUEscience™ is the

world's first laboratory Bluetooth® pH meter and QR buffer system that harnesses the accuracy of standard laboratory grade sensors with the computing power of an Android tablet, expanding measurement capabilities, connecting to the cloud and making your readings SMART.

Technological Innovation

The TRUEscience system offers advanced features which would usually only be seen in more expensive benchtop pH meters such as logging, data storage, multiple electrode connection, alarms and user management at an affordable price.

Connected -

Using low energy Bluetooth® 4.0, the electrode communicates results wirelessly to the tablet in real time. This eliminates cables which can limit truly flexible testing in the field or on-site. Up to 6 SMART Caps can communicate with one tablet at a time, making measurement of different tests or site locations easy.

Credible -

The unique QR buffer system means you always have full traceability on calibrations. See a full history of who has calibrated, which buffers were used and what the results were at any time. Get MSDS and Certificates of Anlaysis for all reagents straight from the TRUEscience App for simple health and safety and quailty checks.

Accessible -

The truly open system lets you use any standard electrode with an S7 connector. As well as removing the need to replace the whole unit every time, this means you can also select an electrode which is specialised for your application, like a spear tip or a flat ended electrode.

Whether you're controlling product quality, testing your site discharge for environmental compliance, checking up on your research in the lab or teaching the next generation of Scientists, you can now measure pH the SMART way.

making measurement SMART

Accessible

Measure Anywhere To and with an effective ra meters you are now fr measurements easily f locations than ever be



Logging The SMART Cap can log continuously and independently tablet and can store up to 10,00 readings before needing to dow The LED flashing indicator lets y at a glance whether the unit is a reading or logging.



Open Electrode System TI three TRUEscience pH elect choose from including glas double junction and flat o use any other standard pH with an S7 connector.

WireFree and Sta TRUEpH is easy. Th wires or cumberso stands and you sin Cap to the side of tubes or rest it on th Stand.



Traceability The innovative QR buffer system automatically loads batch info, expiry date and re-order info into the App and even allows you to access MSDS sheets and batch calibration certifcates for each buffer used.



Full Calibration Suite

Not only does the TRUEpH meter

tell you when it was last calibrated,

it keeps a full calibration history by

user and alerts you to recalibrate using its simple step by step process at intervals you control. It will even alert you if

your buffers have expired.

Credible

pH 10.00 @20°C



View your log a graph within and s







Perfectly Conditioned The unique TRUEscience Storage Stand is the perfect resting place for your electrode keeping the membrane ready to use with minimal solution waste.

Product Ordering Information

Order Code



SMART Cap (with TRUEscience combination pH electrode double junction)

SMART Cap with TRUEscience combination double junction pH electrode/ Soaker bottle Lithium battery, Spare battery covers colour coded (Set of 3)

1201690 (White) **1201822** (Grey)

1201689 (White)

1201821 (Grey)

1201575

1201574

1201711

1144960

1145011





SMART Cap (without pH electrode)

SMART Cap (no pH electrode supplied) Lithium battery, Spare battery covers colour coded (Set of 3)



SMART Cap Complete pH meter Kit

SMART Cap with TRUEscience combination double junction pH electrode/ Soaker bottle1201691(White)USB Port charger Stand for TRUEscience tablet, Bluetooth® enabled Android tablet 7"/8"1201823(Grey)Storage Stand, TRUEscience QR buffers 250ml pH 4, 7 & 10, Spare battery covers colour coded x 3, Lithium battery



SMART Cap pH Meter Kit (without pH electrode)

SMART Cap (no pH electrode supplied), USB Port charger Stand for TRUEscience tablet1201692(White)Bluetooth® enabled Android tablet 7"/ 8" Storage Stand, TRUEscience QR buffers 250ml pH 4, 7 & 101201824(Grey)Spare battery covers colour coded x 3, Lithium battery



TRUEscience pH electrode:

Combination pH electrode glass with S7 connector Combination pH electrode double junction flat with S7 connector Combination pH electrode double junction epoxy with S7 connector

Accessories:

| Replacement Lithium battery for SMART Cap Spare battery covers colour coded Grey, Teal, Red (Set of 3) Spare battery covers colour coded White, Teal, Red (Set of 3) Tablet 7"/ 8" Bluetooth® enabled Android Storage Stand (includes storage tubes) Spare storage tubes for storage stand Charging tablet stand includes USB Cable Height adjustable tablet stand, Stainless steel, non-charging Replacement SMART Clip | | | 1201564 1201710 1201813 1201708 1201695 1201696 1201562 1201563 1201698 |
|--|----------------|----------------|---|
| Buffer solution pH 4.00 \pm 0.01 @ 20°C Red | 250ml | 500ml | 1000ml |
| | 1201459 | 1144885 | 1144887 |

Buffer solution pH 7.00 \pm 0.01 @ 20°C Yellow12014391144803Buffer solution pH 7.00 \pm 0.01 @ 20°C Yellow12014581144955Buffer solution pH 10.00 \pm 0.01 @ 20°C Blue12014601145008

order online at: www.truescience.co.uk



TECHNICAL SPECIFICATION

pH Range **Resolution** Meter Accuracy **Calibration points Buffer Recognition** Automatic temperature compensation pH Electrode (where included) **Electrode** accuracy **Electrode diagnostics** mV range mV accuracy **Temperature accuracy** Store values **Battery Life Battery Type** Tablet and App requirements **Dimensions of Cap** Temperature probe length Weight **Operating/Storage Temperature Calibration Certainty Electrode health reminders** Colour keyed battery covers Multiple Cap Measurements Consumables ordering online Any S7 connector pH electrode Interval logging

8

Complete TRUEpH Package

0.01pH , 0.1mV, 0.1°C +/- 0.002 up to 5 points **QR** Code Input 0°C - 100.0°C Epoxy double junction combination pH +/- 0.02pH Yes +/- 1000mV +/- 0.1mV +/- 0.5°C 10,000 measurements/log 1 year normal use 1/2 AA Lithium Thionyl Chloride, 3.6V, 1200mA See below or specified list of tested models on our website L218 x W48 x D24 mm including temperature probe 130mm long 91.5 g 0-55°C Yes, App indicates calibration status, time and date Yes Yes Yes, up to 6 simultaneously Yes Yes Yes

Tablet Requirements:

The TRUEscience App has been optimised for use with 7 and 8 inch screen Android tablets. The tablet must have Bluetooth® 4.0 or later and Android 5.0 or later operating system. You will require a wifi connection for data sharing and full App functionality. The App will run in the background and will require access to GPS, Camera, Bluetooth® and sharing Apps such as GMail, Dropbox and Googledrive.

0-14 pH units



EN61236-1: 2006 - Annex A 2008: 1999/5/EC Annex III Electrostatic discharge may cause the unit to power off

The Bluetooth® transceiver device meets the requirements for modular transmitter approval as detailed in FCC public notice DA00-1407

www.truescience.co.uk

+44(0)1954 233144 info@truescience.co.uk